


## Introduction



In cases of suspected COVID-19, ASE guidelines suggest that ECGs need not be monitored in order to reduce exposure and contamination [1]. However, as some automated methods for cardiac quantification rely on accurate heart rates (HR) being present, it is important to be able to estimate the heart rates in cases where ECG monitoring was not performed. We propose and evaluate several novel heart rate estimation methods that apply directly to echocardiograms of the left ventricle (LV).

## Methods and Results



### Overview

- A sequence of images is reduced to one-dimensional time series using one of three different dimensionality reduction methods
- U-Net: A convolutional neural network is trained to segment echo images through time [2]. The segmentation area is the signal.
- NMF: The images are factorised into non-negative matrices. The differences of these factors is the signal [3].
- PCA: The images are reduced to principal components, and the magnitude of the first component is the signal.
- The one-dimensional time series is passed through one of two period extraction algorithms
- The heart rate is calculated from the period and the time between frames



### Data

778 DICOM image clips with at least two cardiac cycles, collected from stress echo acquisitions at St George's Hospital, UK, are used for testing. DICOMs encompass a large variety of heart rates, stress echo protocols (both exercise and pharmacological) and view acquisitions (A2C, A4C or SAX-PM).


View/Stage	Contrast		Non-contrast	
	Rest	Stress	Rest	Stress
A2C	76	80	61	57
A4C	76	80	58	57
SAX-PM	61	63	57	52

**Table 1: Breakdown of View/Stage/Contrast of image clips used**

Each combination of methods in each step was applied to each image clip to estimate the heart rate. The estimated heart rate is then compared to the heart rate encoded into the meta data of each DICOM. The method fails if a period could not be extracted from the signal.


Failures	Correlation	Coherence
U-Net	0 (0%)	8 (1%)
NMF	60 (8%)	147 (19%)
PCA	51 (7%)	142 (18%)

**Table 2: Number (Percentage) of failures for each combination**



### Results

U-net dimensionality reduction followed by autocorrelation period extraction proved most accurate for estimating heart rate. The results for this combination are shown in Regression and Bland Altman plots below.



Within 20bpm	Correlation	Coherence
U-Net	97% (753)	96% (744)
NMF	87% (680)	93% (721)
PCA	88% (685)	92% (718)

Within 10bpm	Correlation	Coherence
U-Net	91% (707)	90% (703)
NMF	82% (637)	89% (691)
PCA	82% (638)	88% (688)

**Table 3: Percentage (Number) of clips with estimated heart rate within 20 bpm of actual heart rate**

**Table 4: Percentage (Number) of clips with estimated heart rate within 10 bpm of actual heart rate**

### References

- [1] Mitchell C, Collins K, Hua L, et al. Specific Considerations for Sonographers When Performing Echocardiography during the 2019 Novel Coronavirus Outbreak: Supplement to the American Society of Echocardiography Statement. *J Am Soc Echocardiogr.* 2020;33:654-7
- [2] Upton R, Mumith A, Beqiri A, et al. Automated Echocardiographic Detection of Severe Coronary Artery Disease using Artificial Intelligence. *In review.* 2021.
- [3] Yuan B, Chitturi SR, Iyer G, et al. Machine learning for cardiac ultrasound time series data. *Proc SPIE, Medical Imaging.* 2017;10137.
- [4] Lindström J, Kokko H, Ranta E. Detecting Periodicity in Short and Noisy Time Series Data. *Oikos.* 2016;78:406-10

## Conclusions

We have developed a robust and accurate method of estimating heart rates from multiple echocardiographic views in the presence or absence of contrast, which works even when there is motion of the probe during the acquisition.